
Iumina®








limited in quantity and quality [5-7].

Massively parallel sequencing (MPS) is a high throughput technology that can rapidly generate high quality sequence from targeted areas of the human genome. The technology has reached a level of robustness such that it can be considered a viable approach to analyze challenging forensic samples. Illumina previously released two protocols for sequencing mtDNA (whole mt genome and d-loop region) based on Nextera[®] XT DNA Sample Preparation Kit and the MiSeq FGx™ platform, and published findings in King et al. [8].

Now, Illumina is adding analysis workflow applications on Illumina's cloud-computing environment, BaseSpace™, specifically suited to the mtDNA samples of the Forensic Genomics community. After streaming sequencing data to BaseSpace from your MiSeq FGx™, the data is readily available for data analysis and visualization with the

viewed in the mtDNA Variant Analyzer which provides a modern user interface that enables the user. Multiple samples can be inspected in aggregate, or compared to each other with a responsive experience. When the analyst's trained eye disagrees with the computed call for a locus, the call can be changed and saved back to the cloud. An output report summarizes the variants across samples as well as providing more detailed information per sample when necessary.

The mtDNA work presented herein represents an end-to-end solution for the preparation, analysis, visualization, and reporting supporting today's forensic scientists.



Samples

- * Buccal swabs from volunteers
- * Extracted using QIAamp[®] DNA Investigator kit (Qiagen, Hilden, Germany) according to manufacturer's instructions
- * Quantity, size and quality of libraries assessed with Agilent DNA 1000 kit on the Agilent 2100 Bioanalyzer system

🖊 D-Loop Protocol 🛛 –



- * Human mtDNA D-Loop Hypervariable Region Guide [i]
- * Sequenced using a 2 x 151 cycle run with dual index

Whole Genome Protocol

- * Human mtDNA Genome Guide [ii]
- * Sequenced using a 2 x 251 cycle run with dual index

Improved Human Mitochondrial DNA Analysis using **Next-generation Sequencing and Cloud-based Computing and Storage**

Carey Davis, Kevin Rhodes, Nathalie Mouttham, Anthony Rensfield, Aprajita Mathur, Dan Sa, Van Le-Pham, Amanda Young, Joseph Varlaro, John Walsh

Figure 1. High-level user workflow from sample preparation to reporting of mtDNA samples

- * Normalized PCR amplification products prepared using the Nextera[®] XT DNA Sample Preparation kit
- * Samples were indexed to allow for pooling and demultiplexing
- * Sequencing performed on a MiSeq FGx with MiSeq Reagent Kit v3
- * Data streamed to BaseSpace cloud computing platform during sequencing

Amplicon Start	Amplicon End
29	285
172	408
15997	16236
16159	16401
Table 1. Amplicon locations for the D-Loop Protocol	

Amplicon Start	Amplicon End
397	1892
5195	9796
Table 2. Amplicon locations for the Whole Genome Protocol	

Alignment and variant calling was performed using a bioinformatics pipeline developed and optimized for mtDNA samples from the Forensic Genomics community. This workflow is available as the mtDNA Variant Processor app on Illumina's BaseSpace cloud computing platform. The app supports the processing of both the mtDNA D-Loop and the whole genome samples prepared as specified above. Additionally, custom primers can be used if provided during analysis set-up.

Nextera[®] Adapter Trimming

- * Nextera[®] adapters are removed from the forward and reverse reads
- * Trimming is repeated until no more than 3 adapter bases found on the end of the read
- * Reads with excessive trimming or result in very short amplicons are discarded entirely

Alignment

- * Alignment is performed with BWA-MEM [9]
- * Parameters optimized for single nucleotide homopolymeric regions (i.e. C-stretches)
- * Circular alignment is handled across the origin by identifying the true start and end of reads
- * Indels are realigned to improve alignment and shift to a 3' alignment, particularly important in C-stretches

Account for Primers

* The supplied manifest is used to

- identify the primers and amplicons
- * Primer contributions are removed from the reads for accurate variant calling

Variant Calling

For each position, data is piled up to aid in identifying a consensus call for the position

- * Bases are quality score filtered prior to use in calling
- * A minimum read count is used to avoid positions with practically no coverage or noise
- * A score is calculated for each called position that accounts for the proportion of reads below BaseQ, MapQ, and Analysis Treshold
- * Calls and reporting of alleles are strictly threshold based:

- Percentage of reads for a nucleotide exceeds the Analysis Threshold (AT), it is reported in the output

- Percentage of reads for a

* Ambiguous indels are reported when there is significant data

(exceeds the IT) that supports the indel and doesn't support it (See Figure)


Code	D
А	Adenine
С	Cytosine
Т	Thymine
G	Guanine
R	A and G
Y	C and T
S	G and C
W	A and T
К	G and T
Μ	A and C
В	C and G and T
D	A and G and T
Н	A and C and T
V	A and C and G
Ν	A and C and G and ⁻
-	Deletion
	No Call

[4] M. Richards, V. Macaulay, A. Torroni, H.-J. Bandelt. In search of geographical patterns in European mitochondrial DNA. The American Journal of Human Genetics. 71 (2002) 1168-74. [5] P. Gill, P.L. Ivanov, C. Kimpton, R. Piercy, N. Benson, G. Tully, et al. Identification of the remains of the Romanov family by DNA analysis. Nat Genet. 6 (1994) 130-5. [6] M.M. Holland, T.J. Parsons. Mitochondrial DNA sequence analysis-validation and use for forensic casework. Forensic Sci Rev. 11 (1999) 21-50.

[7] M.R. Wilson, J.A. DiZinno, D. Polanskey, J. Replogle, B. Budowle. Validation of mitochondrial DNA sequencing for forensic casework analysis. Int J Leg Med. 108 (1995) 68-74. [8] J. L. King, B. L. LaRue, N. Novroski, M. Stoljarova, S. B. Seo, X. Zeng, et al., High-quality and high-throughput massively parallel sequencing of the human mitochondrial genome using the Illumina MiSeq, Forensic Sci. Int. Genet. 12 (2014) 128-135. [9] Li H. (2013) Aligning sequence reads, clone sequences and assembly contigs with BWA-MEM. arXiv:1303.3997v1 [g-bio.GN].

nucleotide exceeds the Interpretation Threshold (IT), it is used in the call





de codes for unambiguous consensus calls

Figure 3. An example illustrating the presence of ar ambiguous deletion at position 501. There are an equal number of reads supporting a deletion call and a reference call of G so we report an ambiguous call of g as the consensus.

Code	Definition
а	Ambiguous indel with adenine
С	Ambiguous indel with cytosine
t	Ambiguous indel with thymine
g	Ambiguous indel with guanine
r	Ambiguous indel with A and G
У	Ambiguous indel with C and T
S	Ambiguous indel with G and C
W	Ambiguous indel with A and T
k	Ambiguous indel with G and T
m	Ambiguous indel with A and C
b	Ambiguous indel with C and G and T
d	Ambiguous indel with A and G and T
h	Ambiguous indel with A and C and T
V	Ambiguous indel with A and C and G
n	Ambiguous indel with A and C and G and T

Table 4. Codes for ambiguous consensus calls of insertions and deletions

Output


* VCF file formatted with extra data for use with the mtDNA Variant Analyzer

* BAM file for investigation of the alignment with third party tools (e.g. IGV

while providing all relevant information to understand the data. Additionally, the mtDNA Variant Analyzer app enables simultaneous data analysis across runs and projects for the side-by-side comparison of results collected over time or within a run.

Comparison to rCRS

- * Each sample's consensus genome is displayed against the rCRS genome
- * Navigation options allow to skip to just hotspots (e.g. variants, interpretation threshold warnings, modified calls)



LINKS [i] Illumina http://support.illumina.com/downloads/human_mtdna_d_loop_hypervariable_region_guide_15034858.ilmn [ii] http://support.illumina.com/downloads/human_mtdna_genome_guide_15037958.html

* Different view filters enable the restriction of the view to just hotspots or the D-loop

- view

Position level information

- * For each called position, the coverage, score, and nucleotide distribution can be inspected
- * When comparing samples, the reference and selected sample are shown together

Copyright 2016 Illumina, Inc. All Rights Reserved. Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPro, CytoChip, DesignStudio, Epicentre, GAIIx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, ForenSeq, MiSeq, MiSeqDx, MiSeqFGx, NeoPrep, Nextera, NextBio, NextSeq, Powered by Illumina, SeqMonitor, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

for both the reference and selected sample.