illumına[®]

ForenSeq[™] Universal Analysis Software

Interrogate the broadest range of forensically relevant loci with a simple user interface and powerful analysis algorithms.

Highlights

- Complete, Sample-to-Answer Workflow Manage the sequencing workflow from run set-up to data analysis and report generation
- Perform Automated Data Visualization and Reporting with Simple Graphical User Interface
 Compare samples automatically, generate population statistics, and view data through an intuitive user interface
- Generate Investigative Leads Through Analysis
 Software Advances
 Estimate visible traits such as hair color and eye color as well
 as biogeographical ancestry

Introduction

Illumina ForenSeq Universal Analysis Software is part of the MiSeq FGx[™] Forensic Genomics System—a fully validated¹ nextgeneration sequencing (NGS) system specifically designed for use in forensic genomics applications. The MiSeq FGx System and ForenSeq Universal Analysis Software enable the simultaneous analysis of approximately 200 forensically relevant single nucleotide polymorphisms (SNPs) and short tandem repeats (STRs)—including marker sets not routinely available with traditional methods (Table 1). From sequencing library⁻ and file management to complete data analysis and reporting, the analysis software provides the full functionality needed to analyze and interpret the widest range of casework and database samples.

To maximize flexibility and scalability of the MiSeq FGx System for each laboratory, the ForenSeq Universal Analysis Software is sold separately on a standalone, dedicated server.

Complete, Sample-to-Answer Workflow

The MiSeq FGx Forensic Genomics System uses a pipeline of software applications to perform sequencing runs and complete data analysis (Figure 1). ForenSeq Universal Analysis Software guides the sequencing workflow from run-setup and sample data entry to the final stages of data analysis and report generation. MiSeq FGx Control Software (MCS), installed on the MiSeq FGx instrument control computer, captures flow cell images, operates the flow cell stage, and controls reagent delivery and temperature. During the run, Real-Time Analysis (RTA) software, performs image analysis, base

ForenSeq Universal Analysis Software	Run Set-Up
MiSeq Control Software	Sequencing Chemistry
MiSeq Control Software	Cycle-by-Cycle Imaging
Real Time Analysis Software	Image Analysis, Basecalling, and Quality Scoring
ForenSeq Universal Analysis Software	Alignment, Allele Calling, Genotyping, and Reporting

Figure 1: ForenSeq Universal Analysis Software Workflow — The MiSeq FGx Forensic Genomics System uses the MiSeq FGx Control Software, Real-Time Analysis, and ForenSeq Universal Analysis Software to operate the sequencing instrument and perform data analysis.

Table 1: Simultaneous Analysis of Forensically Relevant Loci with ForenSeq Universal Analysis Software

Feature	Markers ^a
Global Autosomal STRs	27
Y-STRs	24
X-STRs	7
Identity SNPs	94
Phenotypic SNPs	22
Biogeographical Ancestry SNPs	56
Total Number of Loci ^b	> 200

 SNP and STR chromosome locations can be found in the ForenSeq DNA Signature Prep Kit User Guide (support.illumina.com/downloads/forenseq-dna-signature-prepguide-15049528.html).

b. Over 200 markers analyzed when running primer set B.

calling, and assigns base-by-base quality scores. ForenSeq Universal Analysis Software then initiates the final stages of analysis, including demultiplexing, sequence alignment, allele calling, genotyping, and reporting.

^{*} A sequencing "library" is a collection of amplified DNA fragments from a single DNA sample.

orenSeq™					user@	
/ RUN OVERVIEW						
						+ NEW RUN
Create a Run						
Name	Description	ı			Application	
Case No 001-15	Casework	:			Forensic Genomics	
D IMPORT SA	MPLES	• ADD N	EW SAMPLES		익 add exist	ING SAMPLES
SELECT FILE						
Sample Name	Project Name	i7 Index	i5 Index	Sample Type	Міх Туре	
Positive	Case No 001-15	R701	¢ A501	Positive Amplifica	tion Control 🗘 B	÷ ×
Neg	Case No 001-15	R701	¢ A502	Negative Amplific	ation Control 💠 🛛 B	*
001-15 Item 1a-ec	Case No 001-15	R701	¢ A503	\$ Sample	\$ B	÷ ×
001-15 Item 1a-sp	Case No 001-15	R701	¢ A504	\$ Sample	\$ B	*
RB	Case No 001-15	R701	¢ A505	Reagent Blank	\$ B	¢ 🗶
32 Samples						SAVE RUN

Figure 2: ForenSeq Universal Analysis Run Overview Screen – The Run Overview screen features easy data entry fields for Run Name and Run Description. Sample information, index combinations, and application options can also be imported with a tab delimited .txt file. Sample types include Known Sample, Forensic Sample, Reagent Blank, Negative Amplification Control, Positive Amplification Control, and Human Sequencing Control.

Perform Automated Data Visualization and Reporting with Simple Graphical User Interface

Illumina ForenSeq Universal Analysis Software provides an easy data analysis workflow and user-friendly viewing and reporting features. The graphical user interface has a simple, intuitive design that can be viewed from any laboratory computer through the Google Chrome web browser. The Run Overview screen features easy run-setup, sample information, and index tracking options (Figure 2). Within the analysis software, sequencing run data can be associated with specific, user-defined projects and results can be viewed by run or by project in the Projects screen.

At run completion, ForenSeq Universal Analysis Software offers data visualization through several types of display screens. The Sample Details screen provides a summarized view of SNP and STR locus intensity scatter plots. Below the intensity charts, individual allele calls are displayed, along with convenient icons displayed when specific quality flags are triggered (Figure 3). To view more indepth reports on specific loci, drill down using the Locus Detail pop-up screens. The Locus Detail screens display intensity charts, the base-by-base DNA target sequence, stutter, and more (Figure 4). The analysis software also offers a full suite of run metrics and sample quality controls that can be evaluated after the run. Quality metrics and sample quality controls that to quickly scan the results and determine whether controls fall within the recommended ranges.

Project Detail reports and Sample Detail reports can be generated automatically within the software. Each report can be easily printed or exported in .xlsx format. Password protected user accounts ensure secure access to the MiSeq FGx sequencing instrument and to the analysis server.

Generate Investigative Leads through Analysis Software Advances

Compared to current capillary electrophoresis-based methods, Illumina ForenSeq Universal Analysis Software provides a number of data analysis advantages. These include, the capacity to interrogate a greater number of markers (Table 1) and the ability to recover the maximum amount of useful genomic information from degraded DNA, low quality DNA, or complex mixtures². Furthermore, all markers, including SNPs, autosomal, X-, and Y-STRs, can be analyzed simultaneously using a single, streamlined workflow.

Beyond analysis of autosomal, X-, and Y-STRs, the software enables the analysis of marker sets not routinely available with traditional capillary electrophoresis methods. These include a dense set of identity-informative SNPs (iiSNPs)^{3,4}, which are informative for source attribution, phenotypic-informative SNPs (piSNPs)⁵, which provide estimates of eye color (blue, intermediate, brown) and hair color (brown, red, black, blond), and biogeographical ancestry-informative SNPs (aiSNPs)⁶. Biogeographical ancestry estimation is presented as



Figure 3: ForenSeq Universal Analysis Software Sample Details Screen — The Sample Details screen displays overall sample results as an intensity plot, along with detailed summary tables for all loci included in the sequencing library.



Figure 4: ForenSeq Universal Analysis Software Locus Details Screen — The pop-up Locus Detail screen shows the allele call, intensity bar graph, and stutter. View the Repeat Sequence column to see the full, base-by-base target sequence, and any intra-STR variation present in the sample.



Figure 5: ForenSeq Universal Analysis Software Estimation of Visible Traits — The ForenSeq Universal Analysis Software enables estimation of visible traits including eye and hair color. Markers for biogeographical ancestry can also be tested and reviewed.

Table 2: ForenSeq Server

Components and Specifications
$4 \times$ Seagate 1 TB Sata hard drives
RAID controller allowing for data redundancy and speed
Intel 2 GHz \times 64 processor with 6 cores/12 threads
Intel socket R server board
32 GB DDR3 RAM
550 W power supply
ULL, FCC, CE certified
Windows Server 2012 R2 Standard with 5 CALs

Ordering Information

Product	Catalog No.
ForenSeq Universal Analysis Software and Server	SE-550-1001

a principal component analysis (PCA) plot relative to major population groups⁷ (Figure 5). Centroids within the plots then provide logical groupings that give meaningful context to the unknown sample.

Additionally, ForenSeq Universal Analysis Software can detect intra-STR sequence variation that discloses additional identity informative data beyond length variation (Figure 4). These data can be important in kinship analyses and mixture analysis, where DNA contributors share equi-length alleles that can only be distinguished by comparing their base-by-base sequences.

Summary

ForenSeq Universal Analysis Software contains comprehensive sample management and analytical capabilities, including sample and index management, application-specific workflows, data visualization at sample and locus levels, quality flags to simplify data analysis, and easily exportable reports.

Learn More

To learn more about ForenSeq Universal Analysis Software, visit: www.illumina.com/informatics/sequencing-microarray-data-analysis/ forenseq-universal-analysis-software.ilmn

To learn about the MiSeq FGx instrument, visit: www.illumina.com/systems/miseq-fgx.ilmn

To learn about the ForenSeq DNA Signature Prep Kit, visit: www.illumina.com/products/forenseq-dna-signature-kit.ilmn

To learn about mitochondrial DNA analysis with the MiSeq FGx System, visit:

www.illumina.com/applications/forensics/missing_persons_and_ unidentified_human_remains.ilmn

References

- The full MiSeq FGx System workflow is validated per the Scientific Working Group on DNA Analysis Methods (SWGDAM) guidelines (www.swgdam. org).
- Illumina (2014) ForenSeq DNA Signature Prep Kit Data Sheet (www.illumina. com/products/forenseq-dna-signature-kit.ilmn).
- Kidd KK, Pakstis AJ, Speed WC, Grigorenko EL, Kajuna SL, et al. (2006) Developing a SNP panel for forensic identification of individuals. Forensic Sci Int 164(1): 20–32.
- Sanchez JJ, Phillips C, Børsting C, Balogh K, Bogus M, et al. (2006) A multiplex assay with 52 single nucleotide polymorphisms for human identification. Electrophoresis 27(9): 1713–1724.
- Walsh S, Liu F, Wollstein A, Kovatsi L, Ralf A, et al. (2013) The HIrisPlex system for simultaneous prediction of hair and eye colour from DNA. Forensic Sci Int Genet 7(1): 98–115.
- Kidd KK, Speed WC, Pakstis AJ, Furtado MR, Fang R, et al. (2013) Progress toward an efficient panel of SNPs for ancestry inference. Forensic Sci Int Genet 10: 23–32.
- Phillips C, Prieto L, Fondevila M, Salas A, Gómez-Tato A, et al. (2009) Ancestry analysis in the 11-M Madrid bomb attack investigation. PLoS One 4(8): e6583.

Illumina • 1.800.809.4566 toll-free (US) • +1.858.202.4566 tel • techsupport@illumina.com • www.illumina.com

For Research Use Only. Not for use in diagnostic procedures.

© 2015 Illumina, Inc. All rights reserved. Illumina, ForenSeq, MiSeq FGx, and the pumpkin orange color are trademarks of Illumina, Inc. and/ or its affiliate(s) in the U.S. and/or other countries. Pub. No. 1470-2014-002 Current as of 23 October 2015

