


# VeriSeq PGS Library Prep Protocol Guide

For Research Use Only. Not for use in diagnostic procedures.

Quantify Unpurified SurePlex Products	1
Tagment Input DNA	3
Amplify Tagmented DNA	4
Clean Up PCR	5
Normalize Libraries	6
Pool Libraries for the MiSeq System	8
Acronyms	9
Technical Assistance	11



ILLUMINA PROPRIETARY

Document # 100000000812 v02  
April 2016

illumina®

This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY.

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2016 Illumina, Inc. All rights reserved.

**Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiniSeq, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq**, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

## Quantify Unpurified SurePlex Products

### Prepare 1/10 Dilutions of SurePlex Sample and Controls

- 1 Vortex each sample and control.
- 2 Centrifuge at  $280 \times g$  for 1 minute.
- 3 In a new PCR plate, add 45  $\mu\text{l}$  molecular-grade water to the required wells.
- 4 Add 5  $\mu\text{l}$  sample or control.
- 5 Vortex to mix.
- 6 Centrifuge at  $280 \times g$  for 1 minute.
- 7 Set aside on wet ice.

### Qubit Method

- 1 Prepare the working solution according to the manufacturer instructions.
- 2 Add 10  $\mu\text{l}$  of each standard to 190  $\mu\text{l}$  of working solution.
- 3 Add 10  $\mu\text{l}$  of the 1/10 diluted SurePlex sample and 190  $\mu\text{l}$  working solution to each assay tube. Briefly vortex to mix.
- 4 Incubate the assay tubes for 2 minutes.
- 5 Calculate the concentration of each 1/10 diluted SurePlex sample. Convert the units to  $\text{ng}/\mu\text{l}$ .
- 6 Compare the measured concentration of the samples to the values in the VeriSeq PGS-MiSeq QC Assessment Guide.

### Quant-iT Method

- 1 Prepare the working solution according to the manufacturer instructions.
- 2 Add 190  $\mu\text{l}$  working solution to the microplate wells.
- 3 Add 10  $\mu\text{l}$  of each 1/10 diluted SurePlex sample to separate wells. Create duplicates or triplicates. Pipette to mix.
- 4 Add 10  $\mu\text{l}$  of each  $\lambda$  DNA standard to separate wells. Create duplicates or triplicates. Pipette to mix.
- 5 Add 10  $\mu\text{l}$  of each  $\lambda$  DNA standard to separate wells. Create duplicates or triplicates. Pipette to mix.
- 6 Measure the fluorescence.
- 7 Calculate the concentration of each 1/10 diluted SurePlex sample. Convert the units to  $\text{ng}/\mu\text{l}$ .
- 8 Compare the measured concentration of the samples to the values in the VeriSeq PGS-MiSeq QC Assessment Guide.

### Template Dilution to 0.2 $\text{ng}/\mu\text{l}$

- 1 Using BlueFuse Workflow Manager, enter the calculated dsDNA concentration (ng/μl) of the 1/10 diluted SurePlex sample concentration into the VeriSeq PGS-MiSeq Assay Plate.
- 2 According to the BlueFuse Workflow Manager calculations, add molecular-grade water to a new PCR plate.
- 3 Add 5 μl of the 1/10 diluted SurePlex sample to each well of the plate containing molecular-grade water.
- 4 Vortex, and then centrifuge the plate at  $280 \times g$  for 1 minute.
- 5 Set aside on wet ice.

# Tagment Input DNA

## Procedure

- 1 Label a new PCR plate VTA.
- 2 Calculate the total volume of TD. Divide the volume equally among the wells of a PCR 8-tube strip.
- 3 Add 10  $\mu\text{l}$  TD Buffer to each well.
- 4 Add 5  $\mu\text{l}$  ATM to the wells containing TD Buffer.
- 5 Add 5  $\mu\text{l}$  SurePlex amplification product (diluted at 0.2 ng/ $\mu\text{l}$ ) to each sample well.
- 6 Mix at 1,800 rpm for 1 minute.
- 7 Centrifuge at  $280 \times g$  for 1 minute.
- 8 Place on a thermal cycler with a heated lid and run the program:
  - ▶ 55°C for 5 minutes
  - ▶ Hold at 10°C
- 9 Calculate the total volume of NT buffer. Divide the volume equally among the wells of a PCR 8-tube strip.
- 10 Add 5  $\mu\text{l}$  NT Buffer to each well.
- 11 Mix at 1800 rpm for 1 minute.
- 12 Centrifuge at  $280 \times g$  for 1 minute.
- 13 Incubate at room temperature for 5 minutes.

# Amplify Tagmented DNA

## Preparation

- 1 Save the following program on a thermal cycler with a heated lid:
  - ▶ 72°C for 3 minutes
  - ▶ 95°C for 30 seconds
  - ▶ 12 cycles of:
    - ▶ 95°C for 10 seconds
    - ▶ 55°C for 30 seconds
    - ▶ 72°C for 30 seconds
  - ▶ 72°C for 5 minutes
  - ▶ Hold at 4°C

## Procedure

- 1 Print the sample assay plate layout using the BlueFuse Workflow Manager.
- 2 Arrange the index primers in the TruSeq Index Plate Fixture, as follows:
  - ▶ N701–N712 in columns 1–12
  - ▶ S503 in row A, S504 in row C
- 3 Place the plate on the TruSeq Index Plate Fixture.
- 4 Add index adapters. Change tips between each well.
  - ▶ Add 5  $\mu$ l of each Index 1 (i7) adapter to each column.
  - ▶ Add 5  $\mu$ l of each Index 2 (i5) adapter to each row.
- 5 Add 15  $\mu$ l NPM to each well.
- 6 Mix at 1800 rpm for 1 minute.
- 7 Centrifuge at  $280 \times g$  for 1 minute.
- 8 Place on the thermal cycler and run the saved program.

### SAFE STOPPING POINT

If you are stopping, store the sealed plate at -25°C to -15°C for up to 7 days.

# Clean Up PCR

## Preparation

- 1 Thoroughly vortex the AMPure XP beads.
- 2 Prepare fresh 80% EtOH from absolute ethyl alcohol.

## Procedure

- 1 Centrifuge the VTA plate at  $280 \times g$  for 1 minute.
- 2 Add an appropriate volume of beads to a trough.
- 3 Add 45  $\mu\text{l}$  AMPure XP beads to each required well of a clean deep well plate.
- 4 Transfer 45  $\mu\text{l}$  PCR product from the VTA plate to the plate containing beads.
- 5 Mix at 1800 rpm for 1 minute.
- 6 Incubate at room temperature for 5 minutes.
- 7 Pulse centrifuge. *To prevent magnetic bead aggregation, do not centrifuge longer than a pulse.*
- 8 Place on a magnetic stand until the liquid is clear. Keep the plate on the stand during the following steps.
- 9 Discard all supernatant from each well.
- 10 Wash 2 times with 200  $\mu\text{l}$  80% EtOH.
- 11 Remove residual EtOH.
- 12 Air-dry on the magnetic stand for 15 minutes, or until beads are completely dry.
- 13 Add 50  $\mu\text{l}$  RSB to each well.
- 14 Remove the plate from the magnetic stand.
- 15 Mix at 1800 rpm for 1 minute.
- 16 Centrifuge at  $280 \times g$  for 1 minute.
- 17 Place on a magnetic stand until the liquid is clear.
- 18 Transfer 45  $\mu\text{l}$  of each supernatant from each well to a new PCR plate.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  for up to 7 days.

# Normalize Libraries



## WARNING

This set of reagents contains formamide, an aliphatic amide that is a probable reproductive toxin. Personal injury can occur through inhalation, ingestion, skin contact, and eye contact. Dispose of containers and any unused contents in accordance with the governmental safety standards for your region. For more information, see the SDS for this kit, at [support.illumina.com/sds.ilmn](http://support.illumina.com/sds.ilmn).

## Preparation

- 1 Prepare 0.1 N NaOH, as follows.
  - a Add 9.9 ml molecular-grade water to a 15 ml conical tube.
  - b Add 0.1 ml of 10 N NaOH (1/100 dilution).
  - c Vortex to mix.

## Procedure

- 1 Prepare the LNA1/LNB1 mix according to the number of reactions

Reagent	1 rxn	12 rxn*	24 rxn*	48 rxn*	96 rxn*
LNA1	36.8 µl	550 µl	1100 µl	2200 µl	4400 µl
LNB1	8.2 µl	100 µl	200 µl	400 µl	800 µl

\*Includes 20% excess.

- 2 Thoroughly vortex until LNA1/LNB1 mix is homogenized.
- 3 Label a new deep well plate LNP.
- 4 Pour the LNA1/LNB1 mix into a trough.
- 5 Transfer 45 µl LNA1/LNB1 mix to each well.
- 6 Add 20 µl dsDNA to each well.
- 7 Mix at 1800 rpm for 30 minutes.
- 8 Pulse centrifuge. *To prevent magnetic bead aggregation, do not centrifuge longer than a pulse.*
- 9 Place on a magnetic stand until the liquid is clear. Keep the plate on the stand during the following steps.
- 10 Remove and discard all supernatant from each well.
- 11 Discard the tips in an appropriate hazardous waste container. Change tips between samples.
- 12 Wash 2 times, as follows.
  - a Add 45 µl LNW1.
  - b Remove from the magnetic stand.
  - c Shake at 1800 rpm for 5 minutes.
  - d Briefly centrifuge at 280 × g.
  - e Place on a magnetic stand until the liquid is clear.
  - f Remove and discard all supernatant.
- 13 Add 30 µl 0.1 N NaOH to each well.



- 14 Remove from the magnetic stand.
- 15 Mix at 1800 rpm for 5 minutes.
- 16 Make sure that the contents of each well is resuspended.
- 17 Centrifuge at  $280 \times g$  for 1 minute.
- 18 Place on a magnetic stand until the liquid is clear.
- 19 Add 25  $\mu$ l of LNS1 to each well of a new PCR plate.
- 20 Transfer 25  $\mu$ l of supernatant from the LNP plate to the new PCR plate containing LNS1.
- 21 Vortex, and then centrifuge at  $280 \times g$  for 1 minute.

### **SAFE STOPPING POINT**

If you are stopping, store the sealed plate at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  for up to 7 days.

# Pool Libraries for the MiSeq System

## Preparation

- 1 Make sure that the MiSeq System is ready for use.
- 2 If needed, thaw the plate at room temperature, and then vortex.
- 3 Save the following POOL program on a thermal cycler with a heated lid:
  - ▶ 96°C for 3 minutes
  - ▶ 4°C for 5 minutes
  - ▶ Hold at 4°C

## Procedure

- 1 Centrifuge the plate at  $280 \times g$  for 1 minute.
- 2 According to the sample sheet, transfer 5  $\mu\text{l}$  of each normalized library to pool into a LoBind tube.
- 3 Vortex and centrifuge the pooled library.
- 4 Transfer 15  $\mu\text{l}$  library pool to a new PCR tube or PCR 8-tube strip.
- 5 Add 85  $\mu\text{l}$  HT1.
- 6 Record the volumes of library pool and HT1 dispensed in the table below.

Reagent	Recommended Volume	Actual Volume
Pool library	10 $\mu\text{l}$	
HT1	90 $\mu\text{l}$	
Total	100 $\mu\text{l}$	100 $\mu\text{l}$

- 7 Gently vortex and centrifuge the pool/HT1 mixture.
- 8 Place on the preprogrammed thermal cycler and run the POOL program.
- 9 Transfer 600  $\mu\text{l}$  of HT1 into a second clean LoBind tube. Set aside in an ice-water bath.
- 10 Transfer 100  $\mu\text{l}$  of denatured pool/HT1 mixture to the LoBind tube with HT1. Set aside on wet ice.
- 11 Sequence your library according to the *MiSeq System Guide* (document # 15027617).

## Acronyms

Acronym	Definition
ATM	Amplicon Tagment Mix
HT1	Hybridization Buffer
LNA1	Library Normalization Additives 1
LNB1	Library Normalization Beads 1
LNS1	Library Normalization Storage Buffer 1
LNP	Library Normalization Plate
LNW1	Library Normalization Wash 1
NPM	Nextera PCR Master Mix
NT	Neutralize Tagment Buffer
PGS	Pre-Implantation Genetic Screening
RSB	Resuspension Buffer
SCT	Single Cell Tagment
SLB	Single Lysis Buffer
TD	Tagment DNA Buffer
VTA	VeriSeq Tagment Amplicon Plate

## Notes

## Technical Assistance

For technical assistance, contact Illumina Technical Support.

**Table 1** Illumina General Contact Information

<b>Website</b>	www.illumina.com
<b>Email</b>	techsupport@illumina.com

**Table 2** Illumina Customer Support Telephone Numbers

Region	Contact Number	Region	Contact Number
North America	1.800.809.4566	Japan	0800.111.5011
Australia	1.800.775.688	Netherlands	0800.0223859
Austria	0800.296575	New Zealand	0800.451.650
Belgium	0800.81102	Norway	800.16836
China	400.635.9898	Singapore	1.800.579.2745
Denmark	80882346	Spain	900.812168
Finland	0800.918363	Sweden	020790181
France	0800.911850	Switzerland	0800.563118
Germany	0800.180.8994	Taiwan	00806651752
Hong Kong	800960230	United Kingdom	0800.917.0041
Ireland	1.800.812949	Other countries	+44.1799.534000
Italy	800.874909		

**Safety data sheets (SDSs)**—Available on the Illumina website at [support.illumina.com/sds.html](http://support.illumina.com/sds.html).

**Product documentation**—Available for download in PDF from the Illumina website. Go to [support.illumina.com](http://support.illumina.com), select a product, then select **Documentation & Literature**.



Illumina

5200 Illumina Way

San Diego, California 92122 U.S.A.

+1.800.809.ILMN (4566)

+1.858.202.4566 (outside North America)

[techsupport@illumina.com](mailto:techsupport@illumina.com)

[www.illumina.com](http://www.illumina.com)