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<https://www.genomicsengland.co.uk/>

Population Sequencing

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- Conducting the 100,000 Genomes Project
- Sequencing genomes from NHS subjects, with a focus on rare disease and common cancers

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- Re-process 5,000 WGS from individuals with rare disease using latest reference genome (GRCh38)

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Process (align and call variants) samples with latest reference using:

- DRAGEN Germline Pipeline
- DRAGEN Joint Genotyping Pipeline

“Genomics England has found working with Illumina and formerly Edico, in the development and implementation of DRAGEN for our realignment and variant calling programme, very collaborative and constructive. DRAGEN is the first platform to utilise FPGA cards to optimise alignment, which has accelerated our project significantly, leading to more benefits... Genomics England is finalising work to reprocess all 100,000 genomes using DRAGEN and rolling it out prospectively for our production pipeline.”

-Dr. Liz Edwards, Head of Bioinformatics Partnerships

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High-volume genomics services laboratory

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- Large genomic services lab in Alabama, providing sequencing services to research, academic, and commercial customers globally
- Scaling up sequencing operations

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- With HiSeq™ X Ten, HAIB's output tripled, creating a need to analyze the data rapidly to avoid bottlenecks

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DRAGEN on-premise enabled HAIB to phase out local analysis servers and speed up sample analysis to increase sample throughput and eliminate bioinformatics bottlenecks

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- DRAGEN Germline Pipeline
- DRAGEN Joint Genotyping Pipeline
- DRAGEN CNV Pipeline
- DRAGEN Somatic Pipeline
- DRAGEN RNA Pipeline
- DRAGEN VLRD Pipeline

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- Bio IT World, Children's Hospital Of Philadelphia, Edico Set World Record For Secondary Analysis Speed, October 23, 2017, www.bio-it-world.com/2017/10/23/childrens-hospital-of-philadelphia-edico-set-world-record-for-secondary-analysis-speed.aspx, Accessed September 19, 2018.
- The San Diego Union Tribune, Rady Children's Institute sets Guinness world record, February 12, 2018, www.sandiegouniontribune.com/news/health/rd-no-rady-record-20180209-story.html, Accessed September 19, 2018.
- PrecisionFDA Hidden Treasures Warm Up, <https://precision.fda.gov/challenges/1/view/results>, Accessed September 14, 2018.

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Ultra-Rapid Whole Genome Diagnosis for Critically Ill Newborns in 19.5 Hours

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Guinness World Records® designation for Fastest Genetic Diagnosis for Genomic Medicine

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For Research Use Only. Not for use in diagnostic procedures.

illumina®

Children's Hospital of Philadelphiaの例

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1000_WHG

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Guinness World Records® designation for Fastest Analysis of 1,000 Genomes
Children's Hospital of Philadelphia (CHOP)が最速1000ゲノム解析でギネス認定

For Research Use Only. Not for use in diagnostic procedures.

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Figure 3: Flexibility of DRAGEN pipelines—Each DRAGEN pipeline contains a unique set of steps in accordance with its function. Demonstrated by the DRAGEN Gemma Pipeline above, DRAGEN provides the flexibility to insert a variety of input files and produce a range of output documents, enabling users to customize their experience and produce their desired file format.

Custom References

The DRAGEN Reference Builder, also referred to as a hash table, can be used by the customer to generate a non-human or non-standard reference. Created references can be used as an input

Push button: DRAGEN on BaseSpace Sequence Hub makes it easy for labs of varying degrees of informatics expertise to perform secondary analysis in-house at a low cost.

Command line: DRAGEN on-premise offers a command line interface which can be used for single-command launch with an easy-to-learn Linux based command line interface (CLI) or advanced command line.

BaseSpace Sequence Hub	On-Premise
Push button Simple graphical user interface (GUI) Managed service HIPAA compliance* Workgroup capabilities Easy data sharing	Single command launch Easy-to-learn Linux based command line interface (CLI) Simple command line execution
	Advanced command line Noted back to back jobs New default configurations Best for different applications

Figure 4: Options for DRAGEN Implementation—In BaseSpace Sequence Hub, users can simply select the app, input info and start a run. DRAGEN on-premise uses a command line interface. For novice users, an easy-to-learn and operate command line interface can be used. For more advanced users, an advanced command line interface allows for added customization. *HIPAA compliance applies in the US only with BaseSpace Enterprise.

Scalability

The DRAGEN Platform enables labs to scale operations while keeping costs and turnaround times low. DRAGEN can facilitate the expansion of research capabilities in several ways:

- Keeping up with the NovaSeq™ 8000 System:** A single DRAGEN server can keep up with the output of two NovaSeq 8000 instruments with dual S4 flowcells at full capacity.
- Burst capacity:** During times of high capacity with increased sample volumes, labs can scale to DRAGEN on BaseSpace Sequence Hub for burst capacity. The parallel suite of DRAGEN pipelines makes it possible to transfer analysis into BaseSpace Sequence Hub.
- Expanding operations:** A single DRAGEN Platform can be used to run all DRAGEN pipelines and supported sample types. The speed, accuracy and cost efficiency of DRAGEN enable users to scale up operations without compromising turnaround times or quality of results.
- Exome to genomes:** Ramping from whole-exome sequencing (WES) to whole-genome sequencing (WGS) involves a large increase in generated data. DRAGEN enables customers to easily scale from exomes to genomes without large investments in

As we continue to unlock the power of the genome with new and advanced applications, the amount of data generated from next-generation sequencing (NGS) rapidly expands. In 2018, more than 100 petabytes of data were generated by Illumina systems.

To keep up with the vast amount of data, customers require data analysis tools that can efficiently translate the raw sequencing data into meaningful results without compromising accuracy or cost. Furthermore, to harness the benefits of NGS, organizations that are new to this technology will require easy-to-use solutions that reduce the financial and expertise barriers to adoption.

The Illumina DRAGEN Bio-IT Platform is engineered with tight customer collaboration to address the key pain points associated with analysis of NGS data, developing a highly accurate, ultra-rapid secondary analysis solution that meets the needs of both small research labs and population-scale genomic projects.

About the DRAGEN Platform

The Illumina DRAGEN (Dynamic Read Analysis for GENomics) Bio-IT Platform provides secondary analysis of NGS data from genomes, exomes, and transcriptomes. Fundamental features of the DRAGEN Platform address common challenges in genomic analysis, such as lengthy compute times and massive volumes of data. Without compromising accuracy, the DRAGEN Platform delivers quickness, flexibility, and cost efficiency, enabling labs of all sizes and disciplines to do more with their genomic data.

The DRAGEN Platform is a combined hardware and software solution offering a variety of hardware-accelerated secondary analysis pipelines. DRAGEN's suite of analysis pipelines are engineered to run on field-programmable gate array technology (FPGAs), offering hardware-accelerated implementations of genomic analysis algorithms, including BCL conversion, mapping and alignment, sorting, duplicate marking and haplotype variant calling.

improvements.

The DRAGEN Platform is available on-premise and in the cloud through Illumina's BaseSpace Sequence Hub, and all DRAGEN pipelines can be version controlled.

Accurate Results

DRAGEN Platform implementations are built upon world-class algorithms, and stay up to date to adhere to evolving industry standards and best practices. Exceptional analytical sensitivity and specificity are provided for genomic applications and workflows.

The DRAGEN Platform is engineered to remove biases and other sources of error, ensuring that accuracy is maintained across the board. In the 2017 PrecisionFDA Hidden Treasures – Warm Up Challenge, the DRAGEN Platform received the highest score in five out of six accuracy measures for whole-genome variant calling among platforms that recognized all 50 variants.¹ Improved algorithms in each new version of the DRAGEN Platform are designed to increase accuracy (Figure 1).


Figure 1: High accuracy in single nucleotide variant detection—For variant calling comparison with a popular variant calling platform, reference genomic DNA was sequenced using three Illumina library preparation kits and separately analyzed through the DRAGEN Platform or BWA+GATK. Resulting variant calls were compared to a reference genome truth set to determine which variants were falsely detected false positives (FP) or not detected false negatives (FN).

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<http://bit.ly/draBSblog>

BaseSpace Informatics Blog
by Illumina

Home About Subscribe to RSS

May 29, 2019
by [scatreux](#)
in [DRAGEN Features, Releases](#)
[Leave a comment](#)

Somatic Pipeline Improvements with DRAGEN v3.3

by Severine Catreux – Associate Director, Bioinformatics FPGA Development

Significant accuracy gains and speed improvements with DRAGEN v3.3, released April 2019

The DRAGEN engineering and bioinformatics team is excited to announce a new DRAGEN release, v3.3. The second of several releases scheduled for 2019, DRAGEN v3.3 contains improvements across the many pipeline offerings now supported by the DRAGEN platform. This includes accuracy improvements in the germline and somatic pipelines, new features (e.g. CNV DeNovo calling and RNA quantification) and speed gains (Somatic T/N, BCL conversion).

Please see [DRAGEN v3.3 Release Notes](#) for more details. This blog highlights the significant updates to the DRAGEN Somatic Pipeline for small variants, that are part of the v3.3 release.

As one of DRAGEN's core pipelines, the DRAGEN Somatic Pipeline for small variants is utilized by cancer research institutes around the globe. Expanding on the existing functionality, accuracy and speed of the DRAGEN Somatic Pipeline, the v3.3 release placed a high focus on the somatic tumor/normal WGS mode, producing step-function improvements in both accuracy and speed.

Accuracy Improvements:

During the development cycle for v3.3, the DRAGEN engineering and bioinformatics teams took a deep dive into the DRAGEN Somatic Pipeline tumor/normal mode, strengthening the existing algorithm for accuracy improvements. Specific improvements were made in the genotyping module,

BASESPACE SUITE

- BaseSpace Clarity LIMS
- BaseSpace Sequence Hub
- BaseSpace Variant Interpreter
- BaseSpace Cohort Analyzer
- BaseSpace Correlation Engine

RECENT POSTS

- [Somatic Pipeline Improvements with DRAGEN v3.3](#)
- [Enabling Cancer Interpretation At Scale For The Genomics England 100K Genomics Project](#)
- [Doing more with DRAGEN™ v3.2.8](#)
- [New Sequence Quality Metrics in BaseSpace™ Sequence Hub](#)
- [Putting Your Privacy First](#)

ARCHIVES

Select Month

CATEGORIES

tumor and normal control, to allow for differences in library preparation between the tumor sample and the control sample.

These changes are precursors to further accuracy improvements planned for the DRAGEN v3.4 release, specifically in the area of liquid tumor support, where tumor-in-normal contamination will be taken into account.

Accuracy gains of DRAGEN 3.3 over previous DRAGEN versions (3.2) as well as other pipelines (GATK4 MuTect2 and Strelka2) are shown in the plot below. Gains are measured for both SNVs and indels on most datasets.


Figure 1: Comparison of False-Positives (FP) and False-Negatives (FN) between GATK4, Strelka2, DRAGEN 3.2 and DRAGEN 3.3. Lower values are better.

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A New Era of Informatics with Illumina

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© 2019 Illumina, Inc. All rights reserved.