# Infinium™ Japanese Screening Array-24 v1.0 Consortium BeadChip

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[jp.illumina.com](http://jp.illumina.com)

 [www.facebook.com/illumina](https://www.facebook.com/illumina)

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1. ClinVar Database. [www.ncbi.nlm.nih.gov/clinvar](http://www.ncbi.nlm.nih.gov/clinvar). Accessed August 20, 2020.
2. PharmGKB, The Pharmacogenomics Knowledgebase. [www.pharmgkb.org](http://www.pharmgkb.org). Accessed August 20, 2020.
3. National Human Genome Research Institute. [www.genome.gov/](http://www.genome.gov/). Accessed August 20, 2020.
4. Exome Aggregation Consortium (ExAC) Browser. [exac.broadinstitute.org](http://exac.broadinstitute.org). Accessed August 20, 2020.
5. ACMG Recommendations for Reporting of Incidental Findings in Clinical Exome and Genome Sequencing. [www.ncbi.nlm.nih.gov/clinvar/docs/acmg/](http://www.ncbi.nlm.nih.gov/clinvar/docs/acmg/). Accessed August 20, 2020.
6. University of California, Santa Cruz (UCSC) Genome Browser. [genome.ucsc.edu](http://genome.ucsc.edu). Accessed August 20, 2020.
7. NCBI Reference Sequence Blood Group Antigen Gene Mutation Database. [www.ncbi.nlm.nih.gov/projects/gv/rbc/xslcgi.fcgi?cmd=bgmutsystems](http://www.ncbi.nlm.nih.gov/projects/gv/rbc/xslcgi.fcgi?cmd=bgmutsystems). Accessed August 20, 2020.
8. Catalog of somatic mutations in cancer. [cancer.sanger.uk/cosmic](http://cancer.sanger.uk/cosmic). Accessed August 20, 2020.
9. CPIC. [cpicpgx.org](http://cpicpgx.org). Accessed August 20, 2020.
10. Gene Ontology Consortium. [www.geneontology.org](http://www.geneontology.org). Accessed August 20, 2020.
11. Database of Genomic Variants. [dgv.tcag.ca/dgv/app/home](http://dgv.tcag.ca/dgv/app/home). Accessed August 20, 2020.
12. NCBI eQTL Database. [www.ncbi.nlm.nih.gov/projects/gap/eqtl/index.cgi](http://www.ncbi.nlm.nih.gov/projects/gap/eqtl/index.cgi). Accessed August 20, 2020.
13. The Allele Frequency Database. [alfred.med.yale.edu/alfred/snpSets.asp](http://alfred.med.yale.edu/alfred/snpSets.asp). Accessed August 20, 2020.
14. gnomAD, Genome Aggregation Database. [gnomad.broadinstitute.org](http://gnomad.broadinstitute.org). Accessed August 20, 2020.
15. de Bakker PIW, McVean G, Sabeti PC, et al. A high-resolution HLA and SNP haplotype map for disease association studies in the extended human MHC. *Nat Genet.* 2006;38:1166-1172.
16. Neanderthal Genome Browser. [neanderthal.ensemblgenomes.org/index.html](http://neanderthal.ensemblgenomes.org/index.html). Accessed August 20, 2020.
17. PharmGKB, Clinical Annotation Levels of Evidence. [www.pharmgkb.org/page/clinAnnLevels](http://www.pharmgkb.org/page/clinAnnLevels). Accessed August 20, 2020.
18. RefSeq - NCBI Reference Sequence Database. [www.ncbi.nlm.nih.gov/refseq](http://www.ncbi.nlm.nih.gov/refseq). Accessed May 20, 2020.
19. NCBI Genome Reference Consortium. Version GRCh37. [www.ncbi.nlm.nih.gov/grc/human](http://www.ncbi.nlm.nih.gov/grc/human). Accessed May 20, 2020.
20. Illumina (2012) Infinium FFPE QC and DNA Restoration Kit. Accessed May 20, 2020.

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© 2020 Illumina, Inc. All rights reserved.