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June 27, 2018

© 2015 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the US and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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DesignStudio概要

DesignStudioとは

Design a custom solution using our interactive tool

Start A New Design

What's new in DesignStudio?

Introducing AmpliSeq for Illumina

DesignStudio has been improved to include new AmpliSeq for Illumina designers to walk you through the steps of personalizing your custom amplicon panels. Designers for the following application types are supported:

- AmpliSeq Gene Design
- AmpliSeq Hotspot Design

[Start a new design](#)

To learn more, please read the [release notes](#)

Learn about the available assays

- [AmpliSeq Custom DNA](#)
- [Nextera Rapid Capture Custom](#)
- [TruSeq RNA Expression](#)

Your Recent Designs

ID	NAME	ASSAY TECHNOLOGY	STATUS
145981	Cystic Fibrosis	AmpliSeq DNA Gene	Created
144161	Cancer Hotspot	AmpliSeq DNA HotSpot	Created

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
RNA

Coming Soon!

AmpliSeq for Illumina Custom RNA Panel

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

START DESIGN NEW DESIGNS VIEW REPORTS

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Assay Type
DNA RNA

Assay Technology
Compare Technologies
AmpliSeq for Illumina Gene AmpliSeq for Illumina Hotspot Nextera Rapid Capture

Species
Human (hg19) Human (grch38.p2) Mouse (mm10)

Extended Species
Cow (bosTau7) Chicken (galGal4) Dog (canFam3) Chinese Hamster (criGri1)
Pig (susScr3) Sheep (oviAri3) Corn (AGPv3) Rice (IRGSP-1.0)
Soybean (Glyma1.1) Tomato (SL2.40)

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The stringency setting in **DesignStudio** controls how well amplicons included in the final design solution match their target sequence.

Higher stringency settings yield primer designs that result in highly specific amplicons. The number of on-target sequencing reads is maximized (high specificity) and off-target reads are minimized. High stringency designs may come at the expense of lower target coverage *in silico* and can result in design gaps

Lower stringency settings yield primer designs that result in less specific amplicons. The *in silico* coverage of targets is maximized but low stringency designs may come at the expense of a lower number of on-target reads (low specificity) by generating more off-target reads.

The default recommended stringency is **High**.

	On Target	Coverage	
	High Stringency	Medium Stringency	Low Stringency
	 Only high stringency primers	Slight reduction in primer stringency that will increase coverage and may increase off-target	More relaxed primer stringency that will increase coverage with greater off-target risk
	If the <i>in silico</i> coverage is sufficient, this design is the best option	For better coverage of GC-rich regions	For more coverage of GC-rich regions, and tolerance for higher off-target

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Gene Design

Webinar (145990) - Created [Edit](#)

Select Assay 1 Configure Design 2 Manage Targets 3 Submit Design 4 Review Design 5 Get Pricing 6

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Options: CDS Only Exon Only What is CDS?: 5' UTR Coding DNA Sequence (CDS) 3' UTR

Exon Padding: 25 bp (applies only to targets specified by Gene name) [Download Targets](#)

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Gene Name:

Options: CDS Only Exon Only

Exon Padding: (applies only to targets specified by Gene name)

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Webinar (145990) - Targeted [Edit](#)

Select Assay 1 Configure Design 2 Manage Targets 3 **Submit Design 4** Review Design 5 Get Pricing 6

Manage Targets

Design Type AmpliSeq DNA Gene
Species hg19 (Homo sapiens)
Amplicon Size 275
Stringency High

Add Targets

Gene Coordinate File

Submit an Import File
 .csv .bed
Download file templates: .csv .bed [Select a File](#) [Upload File](#)
[Example chromosome formats](#)

Exon Padding 25 bp (applies only to targets specified by Gene name) [Download Targets](#)

Target Type	Name	Chr	Start	Stop	Added
Gene (CDS)	TP53	chr17	7571719	7578811	
06/13/2018					
Gene (Exon)	CDKN2A	chr9	21967750	21975132	
06/13/2018					
Coordinate	EGFR Ex18	chr7	55241589	55241761	
06/13/2018					
Coordinate	EGFR Ex19	chr7	55242382	55242538	

1 - 11 of 11 items

[Delete ALL Targets](#) [Delete Selected Targets](#)

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Submit Design

Important Note: After submission, your design will no longer be editable

Notify me when the design is complete
Email Address:

I agree to these [Terms and Conditions](#)

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Go to DesignStudio | www.illumina.com

illumina® Your AmpliSeq for Illumina design is completed.

DesignStudio

Dear Illumina Customer,

Your **AmpliSeq for Illumina** design, **Webinar**, is completed. To review this design, visit [DesignStudio](#),

[View Design](#)

◆ **Design Summary**

Design ID:	145990
Solution ID:	IAA6791_182
Design Name:	Webinar
Number of Amplicons:	39

Coverage
90.05 %

If you have any questions, contact [Illumina Customer Service](#)

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Modify Design

Modify design will create a new project with settings and targets from this project as the starting point.

Project Name

① Webinar-3

Project Description

New Project Description

Design Stringency

② high [Learn more about stringency](#)


③ Maximize number of design pools for best coverage [?](#)

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Review Design

Solution ID	IAA6791_182
Assay Version	AmpliSeq DNA Gene
Species	Homo sapiens (hg19)
Sample Type	regular
Max Amplicon Length	275
Stringency	high
Number of Amplicons	39
Total target size (bp)	5293
Covered bases (bp)	4766
Number of Pools	2
Amplicons per Pool	21 / 18

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Finalize Design & Get Pricing

Please review your design details below. To get the most accurate pricing for you, we'll need to finalize your design which means that you can no longer alter it.

Design Details

Design Name	Webinar	Species	Homo sapiens (hg19)
Design Id	IAA6791	Amplicons	39
Solution Id	IAA6791_182	Pools	2
Design Type	AmpliSeq DNA Gene	Amplicons per Pool	21 / 18
		Coverage	90.05%

Finalize Design & Get Pricing Cancel

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Webinar (145990) - Finalized

Select Assay Configure Design Manage Targets

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BaseSpace DESIGNSTUDIO SEQUENCING

START DESIGN VIEW DESIGNS VIEW REPORTS ? Help Shinji

View Designs Project Id: [] Customer Email: [] Clear Q Search: []


Legend: Dual Pool

AmpliSeq DNA	Nextera Rapid Capture	TruSeq Custom Amplicon	TruSeq Targeted RNA	All States			Added	Actions
Type	Design ID	Solution ID	Design Name	Amplicon Size	Amplicon Count	Coverage		
Gene	[]	[]	[]	175	1277 Amplicons	99%	04/30/2018	Select
Gene	[]	[]	[]	175	1190 Amplicons	95%	04/30/2018	Select
Gene	[]	[]	[]	175	1255 Amplicons	96%	04/28/2018	Select
Gene	[]	[]	[]	175	1161 Amplicons	92%	04/27/2018	Select
Gene	[]	[]	[]	275	16 Amplicons	99%	04/24/2018	Select
Gene	[]	[]	[]	275	15 Amplicons	94%	04/24/2018	Select
Gene	[]	[]	[]	275	15 Amplicons	94%	04/18/2018	Select
Gene	[]	[]	[]	275	524 Amplicons	99%	03/28/2018	Select
Gene	[]	[]	[]	275	414 Amplicons	98%	Orderable grch38.p2 03/27/2018	Select
Gene	[]	[]	[]	275	632 Amplicons	98%	Orderable grch38.p2 03/27/2018	Select

Select

- Copy
- Transfer
- Delete
- Download
- Manifest**

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* Tomaiuolo, A.C., Sofia, V.M., Surace, C., Majo, F., Genovese, S., Petrocchi, S., Grotta, S., Alghisi, F., Lucidi, V. and Angioni, A. (2015). Relationship between CFTR and CTRC variants and the clinical phenotype in late-onset cystic fibrosis disease with chronic pancreatitis. *J Mol Diagn* 17, 171-178.

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Cystic Fibrosis (146303)- Targeted [Edit](#)

Select Assay Configure Design **Manage Targets** Submit Design Review Design Get Pricing

1 2 **3** 4 5 6

Manage Targets

Design Type AmpliSeq DNA Gene
Species grch38_p2 (Homo sapiens)
Amplicon Size 275
Stringency High

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Add Targets

Gene Coordinate File

Gene Name: [Add Gene](#) For full region design including introns, use Coordinate input.

Options: CDS Only Exon Only What is CDS ? : 5' UTR Coding DNA Sequence (CDS) 3' UTR

Exon Padding: (applies only to targets specified by Gene name) [Download Targets](#)

Target Type	Name	Chr	Start	Stop	Added
Gene (CDS)	CFTR	chr7	117479962	117668664	06/14/2018
Gene (CDS)	SPINK1	chr5	147824579	147839231	06/14/2018
Gene (CDS)	PRSS1	chr7	142749467	142753076	06/14/2018
Gene (CDS)	PRSS2	chr7	142770942	142774564	06/14/2018
Gene (CDS)	CTRC	chr1	15438441	15446658	06/14/2018
Gene (CDS)	CASR	chr3	122183682	122286497	06/14/2018

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COSMIC (Catalogue Of Somatic Mutation In Cancer)

https://cancer.sanger.ac.uk/cosmic

The screenshot shows the COSMIC website interface. At the top, there is a navigation bar with a 'Genome Version' dropdown menu. A red box highlights this menu, with a hand-drawn arrow labeled '1' pointing to it. The dropdown shows 'GRCh37 ✓' and 'GRCh38'. Below the navigation bar, there is a search bar containing the text 'KRAS G12D', highlighted with a red box and a hand-drawn arrow labeled '2'. To the right of the search bar is a red 'SEARCH' button, with a hand-drawn arrow labeled '3' pointing to it. The main content area includes a 'COSMIC v85, released 08-MAY-18' announcement, a 'COSMIC News' section with various updates, and a 'Projects' section listing different COSMIC projects. At the bottom, there is a cookie consent banner with 'Accept' and 'Cookie Preferences' buttons.

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COSMIC
Catalogue Of Somatic Mutations In Cancer

GRCh 37

Projects ▾ Data ▾ Tools ▾ News ▾ Help ▾ About ▾ Genome Version ▾ Search COSMIC... **SEARCH** Login ▾

COSMIC search results

Your search term "**KRAS G12D**" returned results in **2** sections of the database. [More...](#)

Genes (0 hits) **Mutations (3)** SNPs (0) Cancer (0) Tumour Site (0) Samples (0) Pubmed (2) Studies (0)

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COSMIC
Catalogue Of Somatic Mutations In Cancer

GRCh 37

Projects ▾ Data ▾ Tools ▾ News ▾ Help ▾ About ▾ Genome Version ▾ Search COSMIC... **SEARCH** Login ▾

Mutation
COSM521

GRCh37 · COSMIC v85

Overview

This section shows a general overview of the selected mutation. It describes the source of the mutation i.e gene name/sample name/tissue name with unique ID, and also shows the mutation syntax at the amino acid and nucleotide sequence level. You can see more information on our [help pages](#).

Mutation ID COSM521 **ID**

Gene name [KRAS](#)

AA mutation p.G12D (Substitution - Missense, position 12, G→D)

CDS mutation c.35G>A (Substitution, position 35, G→A)

Nucleotides inserted n/a

Genomic coordinates GRCh37, 12:25398284..25398284, view [Ensembl contig](#)

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HomoloGene n/a

Ever confirmed somatic? Yes

FATHMM prediction Pathogenic (score 0.98)

Remark n/a

Recurrent n/a

Drug resistance n/a

[Reset page](#)

- Overview
- Tissue distribution
- Samples
- Pathways affected
- References

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UCSC Genome Browser on Human Feb. 2009 (GRCh37/hg19) Assembly

move <<< << < > >> >>> zoom in 1.5x 3x 10x base zoom out 1.5x 3x 10x 100x

chr7:123,478,452-126,536,741 3,058,290 bp. EGFR go

chr7 (p11.1-q11) 1031.1 1012 1041 1043.44

Scale chr7: 124,000,000 124,500,000 125,000,000 125,500,000 126,000,000 126,500,000 hg19

INMS497_167 Targets
Project 144324 Solution INMS497_167 Targets
Project 144324 Solution INMS497_167 Amplicons
Project 144324 Solution INMS497_167 Gaps

UCSC Genes (RefSeq, GenBank, CCDS, Rfam, tRNAs & Comparative Genomics)
GC Percent in 5-Base Windows

Catalogue of Somatic Mutations in Cancer V82
Repeating Elements by RepeatMasker

UCSC annotations of RefSeq RNAs (mRNAs and ncRNAs)

move start < 2.0 > Click on a feature for details. Click or drag in the base position track to zoom in. Click side bars for track options. Drag side bars or labels up or down to reorder tracks. Drag tracks left or right to new position. Press "?" for keyboard shortcuts. move end < 2.0 >

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UCSC Genome Browser on Human Feb. 2009 (GRCh37/hg19) Assembly

move <<< << < > >> >>> zoom in 1.5x 3x 10x base zoom out 1.5x 3x 10x 100x

chr7:55,177,510-55,275,773 98,264 bp. enter position, gene symbol, HGVS or search terms go

chr7 (p11.2) 21.3 14.3 14.1 21.11 22.1 q31.1 7033 q34 q35

Scale chr7: 55,190,000 55,200,000 55,210,000 55,220,000 55,230,000 55,240,000 55,250,000 55,260,000 55,270,000 hg19

UCSC Genes (RefSeq, GenBank, CCDS, Rfam, tRNAs & Comparative Genomics)
EGFR-AS1
EGFR

EGFR
EGFR
EGFR
EGFR
EGFR
EGFR


GC Percent
70
30
GC Percent in 5-Base Windows

COSMIC Regions
Catalogue of Somatic Mutations in Cancer V82

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UCSC Genome Browser on Human Feb. 2009 (GRCh37/hg19) Assembly

chr7:55,241,612-55,241,738 127 bp.

UCSC Genes (RefSeq, GenBank, CCDS, Rfam, tRNAs & Comparative Genomics)

GO Percent in 5-Base Windows

Drag-and-select

- Hold **Shift+drag** to show this dialog
- Hold **Alt+drag** to add a highlight
- Hold **Ctrl+drag** (Windows) or **Cmd+drag** (Mac) to zoom
- To cancel, press **esc** anytime or drag mouse outside image
- Highlight the current position with **h** then **m**
- Clear all highlights with **View - Clear Highlights** or **h** then **c**

Highlight color: #aaedff [Reset](#)

Don't show this again and always zoom with shift.
Re-enable via 'View - Configure Browser' (c then f)

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