

## Synthesize cDNA

- 1 Vortex RCS1 for 5 seconds.
- 2 Centrifuge at  $600 \times g$  for 5 seconds.
- 3 Dilute according to your input RNA:
  - ▶ Dilute 50 ng intact total RNA with nuclease-free water to 5  $\mu$ l.
  - ▶ Dilute  $\geq 200$  ng degraded RNA with nuclease-free water to 3  $\mu$ l.
- 4 Add diluted RNA to a plate:
  - ▶ Add 5  $\mu$ l diluted intact total RNA to the CDP1 plate.
  - ▶ Add 3  $\mu$ l diluted degraded RNA to the CDP plate.
- 5 Combine the following volumes in a 1.7 ml microcentrifuge tube. Multiply each volume by the number of samples being prepared.

Reagent	Intact Total RN A Volume ( $\mu$ l)	Degraded RN A Volume ( $\mu$ l)
RCS1	4.4	4.4
ProtoScript II Reverse Transcriptase	1.1	2.2
10X DTT (0.1M)*	0	1.1
Total volume per pool	5.5	7.7

\* Included with ProtoScript II Reverse Transcriptase reagent.

- 6 Invert to mix.
- 7 Centrifuge at  $600 \times g$  for 5 seconds.
- 8 Distribute evenly into an 8-tube strip.

- 9 Add the volume for your plate:
  - ▶ Add 5  $\mu$ l to the CDP1 plate.
  - ▶ Add 7  $\mu$ l to the CDP plate.
- 10 Shake at 1600 rpm for 20 seconds.
- 11 Centrifuge at  $280 \times g$  for 1 minute.
- 12 Place on the thermal cycler and run the CDNASYN1 or CDNASYN2 program.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at  $2^{\circ}\text{C}$  to  $8^{\circ}\text{C}$  for up to 7 days. Alternatively, leave on the thermal cycler overnight.

## Hybridize Oligo Pool

- 1 Combine the following volumes in a 1.7 ml microcentrifuge tube. Multiply each volume by the number of reactions being prepared.

Reagent	Volume ( $\mu$ l)
TOP	5.5
Additional TOP or TE buffer	5.5
Total volume per reaction	11

- 2 Vortex for 5 seconds.
- 3 Centrifuge at  $600 \times g$  for 5 seconds.
- 4 Distribute into an 8-tube strip.
- 5 Add 10  $\mu$ l to the CDP or CDP1 plate.
- 6 Shake at 1600 rpm for 20 seconds.
- 7 Incubate at room temperature for 1 minute.
- 8 Vortex OB1 for 5 seconds.
- 9 Add 30  $\mu$ l OB1 to the CDP or CDP1 plate.
- 10 Shake at 1600 rpm for 1 minute.
- 11 Place on the thermal cycler and run the ANNEAL program.
- 12 Centrifuge briefly.

## Wash, Extend, and Ligase Bound Oligos

- 1 Transfer all supernatant to the HYP plate.
- 2 Place on a magnetic stand until liquid is clear.
- 3 Remove and discard all of the supernatant.
- 4 Move from the magnetic stand to a bench.
- 5 Add 100 µl AM1 to each well.
- 6 Shake at 1800 rpm for 2 minutes.
- 7 Centrifuge at 280 × g for 5 seconds.
- 8 Unseal and place on a magnetic stand until liquid is clear.
- 9 Remove and discard all supernatant.
- 10 Move from the magnetic stand to a bench.
- 11 Add 175 µl UB1.
- 12 Shake at 1800 rpm for 2 minutes.
- 13 Centrifuge at 280 × g for 5 seconds.
- 14 Unseal and place on a magnetic stand until liquid is clear.
- 15 Invert ELM4.
- 16 Remove and discard all supernatant.
- 17 Move from the magnetic stand to a bench.
- 18 Add 40 of ELM4.
- 19 Shake at 1800 rpm for 2 minutes.
- 20 Centrifuge at 280 × g for 5 seconds.
- 21 Place on the 37°C preheated microheating system and incubate for 45 minutes.
- 22 Remove the adhesive seal from the plate.
- 23 Unseal and place on a magnetic stand until liquid is clear.
- 24 Remove and discard all supernatant.
- 25 Add 50 µl of UB1.

## Amplify Libraries

- 1 Arrange Index 1 (i7) adapters in columns 1–12.
- 2 Arrange Index 2 (i5) adapters in rows A–H.
- 3 Place the plate on the TruSeq Index Plate Fixture.
- 4 Add 4 µl of each Index 1 adapter down each column.
- 5 Add 4 µl of each Index 2 adapter across each row.
- 6 Remove and discard all supernatant from the HYP plate.
- 7 Remove from the magnetic stand.
- 8 Add 22.5 µl diluted HP3.
- 9 Shake at 1800 rpm for 30 seconds.
- 10 Incubate at room temperature for at least 5 minutes.
- 11 Create the amplification mix:
  - ▶ **96 libraries**—Add 56 µl TDP1 to 2.8 ml of PMM2.
  - ▶ **48 libraries**—Combine 28 µl TDP1 and 1.4 ml PMM2 in a 1.7 ml microcentrifuge tube.
  - ▶ **16 libraries**—Combine 9.2 µl TDP1 and 460 µl PMM2 in a 1.7 ml microcentrifuge tube.
- 12 Invert to mix.
- 13 Add 22 µl to the IAP plate.
- 14 Unseal the HYP plate.
- 15 Place on a magnetic stand until liquid is clear.
- 16 Transfer 20 µl supernatant from the HYP plate to the IAP plate.
- 17 Shake at 1600 rpm for 30 seconds.
- 18 Centrifuge at 280 × g for 1 minute.
- 19 Place on the thermal cycler and run the program.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

## Clean Up Libraries

- 1 Add 85 µl AMPure XP Beads to the CLP plate.
- 2 Centrifuge the IAP plate at 280 × g for 1 minute.
- 3 Unseal the IAP plate.
- 4 Transfer all supernatant to the CLP plate.
- 5 Shake at 1800 rpm for 2 minutes.
- 6 Centrifuge the plate at 280 × g for 5 seconds.
- 7 Incubate room temperature for 15 minutes
- 8 Unseal and place on a magnetic stand until liquid is clear.
- 9 Remove and discard 135 µl supernatant.
- 10 Wash 2 times with 200 µl 80% EtOH.
- 11 Air dry on the magnetic stand for 15 minutes.
- 12 Remove from the magnetic stand.
- 13 Add 15 µl RSB.
- 14 Shake at 1800 rpm for 2 minutes.
- 15 Centrifuge at 280 × g for 5 seconds.
- 16 Return RSB to 2°C to 8°C storage.
- 17 Incubate the plate at room temperature for 2 minutes.
- 18 Unseal and place on a magnetic stand until liquid is clear.
- 19 Transfer 12.5 µl supernatant to the LNP plate.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

## Pool and Quantify Libraries

- 1 Transfer 5 µl from the LNP plate to a 2 ml microcentrifuge tube.
- 2 Vortex for 5 seconds.
- 3 Centrifuge at 600 × g for 5 seconds.
- 4 Load 1 µl pooled library onto the Standard Sensitivity NGS Fragment Analysis Kit or DNA 1000 Kit.
- 5 Determine the concentration of the pooled library.
- 6 Select the **Region Analysis** tab.
- 7 Drag the blue region lines to capture the 100–300 bp region.
- 8 Dilute each pooled library to 4 nM using RSB.
- 9 Denature and dilute the 4 nM library to the concentration for the sequencing instrument you are using. See the denature and dilute guide for your instrument.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C.

## Acronyms

Acronym	Definition
ACD1	Amplicon Control DNA 1
ACP1	Amplicon Control Oligo Pool 1
TSO	TruSight Oligos
CLP	Clean-up Plate
EBT	Elution Buffer with Tris
ELM4	Extension Ligation Mix 4
FPU	Filter Plate Unit
HT1	Hybridization Buffer
HYP	Hybridization Plate
IAP	Index Amplification Plate
LNA1	Library Normalization Additives 1
LNB1	Library Normalization Beads 1
LNP	Library Normalization Plate
LNS2	Library Normalization Storage Buffer 2
LNW1	Library Normalization Wash 1
OHS2	Oligo Hybridization for Sequencing Reagent 2
PAL	Pooled Amplicon Library