

illumina cancer panel portfolio

| Application | | Solid Tumor | | | | | | | |
|-------------------|-----------------|-----------------------------------|--------------------|---|-----------------------------------|--|--------------------------|--|--------------------------|
| Product | | AmpliSeq™ for Illumina BRCA Panel | TruSight™ Tumor 15 | AmpliSeq for Illumina Cancer Hotspot Panel v2 | AmpliSeq for Illumina Focus Panel | AmpliSeq for Illumina Comprehensive Panel v3 | TruSight Tumor 170 | AmpliSeq for Illumina Comprehensive Cancer Panel | TruSight Oncology 500 |
| Genomic Content | # genes | 2 | 15 | 50 | 52 | 161 | 170 | 409 | 523 |
| | Sample Type | DNA | DNA | DNA | DNA, RNA | DNA, RNA | DNA, RNA | DNA | DNA, RNA |
| Input | Sample Input | 20 ng* | 20 ng* | 10 ng* | 10 ng DNA* 10 ng RNA* | 20 ng DNA* 20 ng RNA* | 40 ng DNA* 40 ng RNA* | 40 ng* | 40 ng DNA* 40 ng RNA* |
| | SNV/InDel | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Variants Detected | CNV | | ✓ | | ✓ | ✓ | ✓ | ✓ | ✓ |
| | Fusion | | | | ✓ | ✓ | ✓ | | ✓ |
| | Splice Variants | | | | ✓ | ✓ | ✓ | | ✓ |
| | MSI | | | | | | | | ✓ |
| | TMB | | | | | | | | ✓ |
| | | | | | | | | | |

*FFPE-compatible.

For more information, visit www.illumina.com/cancerpanels

For Research Use Only. Not for use in diagnostic procedures.

Illumina cancer panel portfolio

| Application | | Germline | | Hematology | | Hematology-Solid Tumor | | |
|-------------------|-----------------|----------------------------------|-----------------|------------------|-------------------------------|--|---------------------|-------------------------|
| Product | | AmpliSeq for Illumina BRCA Panel | TruSight Cancer | TruSight Myeloid | AmpliSeq for Illumina Myeloid | AmpliSeq for Illumina Childhood Cancer Panel | TruSight RNA Fusion | TruSight RNA Pan-Cancer |
| Genomic Content | # genes | 2 | 94 | 54 | 69 | 203 | 507 | 1,385 |
| Input | Sample Type | DNA | DNA | DNA | DNA, RNA | DNA, RNA | RNA | RNA |
| | Sample Input | 20 ng | 50 ng | 50 ng | 20 ng DNA 10 ng RNA | 20 ng DNA* 20 ng RNA* | 10–100 ng* | 10–100 ng* |
| Variants Detected | SNV/InDel | ✓ | ✓ | ✓ | ✓ | ✓ | | ✓ |
| | CNV | | ✓ | | ✓ | ✓ | | |
| | Fusion | | | | ✓ | ✓ | ✓ | ✓ |
| | Splice Variants | | | | ✓ | ✓ | | |
| | Expression | | | | | ✓ | | ✓ |
| | | | | | | | | |

*FFPE-compatible

For more information, visit www.illumina.com/cancerpanels

For Research Use Only. Not for use in diagnostic procedures.

©2018 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html. PUB 1170-2014-005 QB 5484

