








1. Click on  to import columns into Data tables (3 choices)

Note: Differential Methylation Profile Table, Sample Methylation Profile Table, or Group Methylation Profile Table (not shown)



2. Select "Column Import"
3. Click "Browse" to select the SNP table file
4. Click "OK" to import SNP information into the active table

After import, three new columns will appear on the left side of the table.

In the example shown, SNP data is imported into the Differential Methylation Profile Table




Figure 5: Importing the Supplemental SNP List into the GenomeStudio Methylation Module —The Supplementary SNP list can be directly imported into the GenomeStudio Methylation Module project tables by following the 4 steps illustrated above.